



6 in 1 Laser Hair Growth Comb Brush



Features :

The 6-in-1 Cosmetic Laser Comb promotes hair growth and metabolism of the scalp by 6 innovative methods:


- Soft laser
- Pulsed Light
- Micro current massage
- Ionic +/- Vibrate
- 12ml space for adding Growth or Treatment Hair Tonic.



Introduction

The 6-in-1 Cosmetic Laser Comb promotes hair growth and metabolism of the scalp by six innovative methods: soft laser, Pulsed Light , Micro current massage, Ionic, Vibrate and 12ml space for adding Growth or Treatment Hair Tonic

Just like any living cell also benefit your hair from light. The most common example of this is photosynthesis, wherein the energy of solar radiation is converted into chemical energy. Similarly, does the laser light therapy. The laser light is radiated directly targeted and painless to the scalp.



Specification:

Voltage: DC 5V

Battery: 3.7V 600mA

Net Weight:166g/piece

Water tank volume:12ml

Vibration frequency: 5000-9000rpm

Contact Info: sales@air-beautyspa.com



Characteristics:

- 1, BIO Micro electric current mainly for fur hair follicle tissue motion.
- 2, Laser can activate the skin cells.
- 3, Portable design for personal home and traveling easy to carry
- 4, Skin Care Ion +/- Cleaning & nourishing restore hair follicles healthy
- 5, Vibration massage
- 6, Support 12ml space for adding Growth or Treatment Hair Tonic

Efficacy:

- Infrared Therapy Activates cells and promotes metabolism, removing dead cells and reproduction new ones.
- Ozone clean up the bacteria in the hair follicle clogging junk, restoring hair follicles healthy. It stimulates the scalp and promotes blood circulation.
- Impulse Acupuncture Adjustment can efficiently activate cells, dilate capillaries and awaken deeply-slept hair follicles.
- With Ionic +/- it can adjust oil secretion, improve hair loss and hair quality
- It can prevent headaches if you often used

Using the unit

Caution:

- Do not use it on broken skin or areas of the skin with rash
- Avoid water on the comb body
- Make sure your hands are dry before you remove or charge combo brush
- Don't open and assembly by yourself




Contact Info: sales@air-beautyspa.com

Using the unit

1. Add hair tonic

- Remove the Water tank button
- Move out the brush head
- Add the tonic at head back side



2. Infrared Therapy


- Red Color: Activates cells and promotes metabolism, removing dead cells and reproduction new ones
- Blue color : Relax and adjust oil secretion
- Mix color: awaken deeply-slept hair follicles.



Contact Info: sales@air-beautyspa.com


3. Ionic Therapy

- Ionic + : Clean dust and dead cell
- Ionic - : Nourishing , you can add growth or treatment hair tonic
- And so on you can switch on the laser and vibration together to make deep massage



4. Micro current massage

- Switch on the vibration mode to make deep massage



Contact Info: sales@air-beautyspa.com